

MEDICAL ELECTRONIC SYSTEMS

WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	090214	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	09_Feb_15	LOT#	011113001	LOT#	090214002	LOT#	090214003		
For use with:	Conc.			Conc.			Conc	MSC	
SQA-V Automated	44	+/-	6.2	22	+/-	4.4	0.0	0.0	
SQA-V Visualization	43	+/-	10.8	23	+/-	5.8	0.0	0.0	
Hemocytometer	44	+/-	11	23	+/-	5.8	0.0	NA	
Makler *	47	+/-	11.8	25	+/-	6.3	0.0	NA	
Fixed Coverslip	45	+/-	11.3	23	+/-	5.8	0.0	NA	

* Makler chamber was loaded with 5ul.