

MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	160315	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	16_March_16	LOT# 160315001			LOT# 160315002			LOT# 160315003	
For use with:		Conc.			Conc.			Conc	MSC
SQA-V Automated		45	+/-	6.3	23	+/-	4.6	0.0	0.0
SQA-V Visualization		45	+/-	11.3	23	+/-	5.8	0.0	0.0
Hemocytometer		45	+/-	11.3	23	+/-	5.8	0.0	NA
Makler		45	+/-	11.3	23	+/-	5.8	0.0	NA
Fixed Coverslip		45	+/-	11.3	23	+/-	5.8	0.0	NA