

MEDICAL ELECTRONIC SYSTEMS

WWW.MES-GLOBAL.COM

| QwikCheck™ Beads Batch Release Form | | | | | | | | | |
|--|-----------|----------------|-----|------|----------------|-----|-----|-------------------------|------------|
| BATCH # | 0111015 | Level 1 | | | Level 2 | | | NEGATIVE CONTROL | |
| EXP DATE: | 11-Oct-16 | LOT# 111015001 | | | LOT# 111015002 | | | LOT# 111015003 | |
| For use with: | | Conc. | | | Conc. | | | Conc | MSC |
| SQA / QwikCheck gold Automated | | 48 | +/- | 6.7 | 25 | +/- | 5.0 | 0.0 | 0.0 |
| SQA-V Visualization | | 49 | +/- | 12.3 | 25 | +/- | 6.3 | 0.0 | 0.0 |
| Hemocytometer | | 48 | +/- | 12.0 | 25 | +/- | 6.3 | 0.0 | NA |
| Makler | | 50 | +/- | 12.5 | 26 | +/- | 6.5 | 0.0 | NA |
| Fixed Coverslip | | 49 | +/- | 12.3 | 25 | +/- | 6.3 | 0.0 | NA |