

MEDICAL ELECTRONIC SYSTEMS

WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	080216	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	8-Feb-17	LOT# 080216001			LOT# 080216002			LOT# 080216003	
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated		49	+/-	6.9	26	+/-	5.2	0.0	0.0
SQA-V Visualization		49	+/-	12.3	28	+/-	7.0	0.0	NA
Hemocytometer		49	+/-	12.3	26	+/-	6.5	0.0	NA
Makler		51	+/-	12.8	26	+/-	6.5	0.0	NA
Fixed Coverslip		49	+/-	12.3	26	+/-	6.5	0.0	NA