

MEDICAL ELECTRONIC SYSTEMS

WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	261216	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	26_Dec_17	LOT# 261216001			LOT# 261216002			LOT# 261216003	
For use with:	Conc.			Conc.			Conc	MSC	
SQA / QwikCheck gold Automated	47	+/-	6.6	24	+/-	4.8	0.0	0.0	
SQA-V Visualization	49	+/-	12.3	26	+/-	6.5	0.0	NA	
Hemocytometer	49	+/-	12.3	26	+/-	6.5	0.0	NA	
Makler	48	+/-	12.0	26	+/-	6.5	0.0	NA	
Fixed Coverslip	48	+/-	12.0	25	+/-	6.3	0.0	NA	