

**MEDICAL ELECTRONIC SYSTEMS**

**WWW.MES-GLOBAL.COM**

<b>QwikCheck™ Beads Batch Release Form</b>									
<b>BATCH #</b>	141117	<b>Level 1</b>			<b>Level 2</b>			<b>NEGATIVE CONTROL</b>	
<b>EXP DATE:</b>	14_Nov_2018	<b>LOT#</b>	141117001	<b>LOT#</b>	141117002	<b>LOT#</b>	141117003		
<b>For use with:</b>	<b>Conc.</b>			<b>Conc.</b>			<b>Conc</b>	<b>MSC</b>	
<b>SQA / QwikCheck gold Automated</b>	48	+/-	6.7	26	+/-	5.2	0.0	0.0	
<b>SQA-V Visualization</b>	47	+/-	11.8	27	+/-	6.8	0.0	NA	
<b>Hemocytometer</b>	50	+/-	12.5	25	+/-	6.3	0.0	NA	
<b>Makler</b>	46	+/-	11.5	24	+/-	6.0	0.0	NA	
<b>Fixed Coverslip</b>	47	+/-	11.8	25	+/-	6.3	0.0	NA	