

MEDICAL ELECTRONIC SYSTEMS**WWW.MES-GLOBAL.COM**

QwikCheck™ Beads Batch Release Form									
BATCH #	210317	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	21_March_2018	LOT#	210317001	LOT#	210317002	LOT#	210317003		
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated		47	+/-	6.6	25	+/-	5.0	0.0	0.0
SQA-V Visualization		46	+/-	11.5	26	+/-	6.5	0.0	NA
Hemocytometer		49	+/-	12.3	25	+/-	6.3	0.0	NA
Makler		47	+/-	11.8	26	+/-	6.5	0.0	NA
Fixed Coverslip		46	+/-	11.5	26	+/-	6.5	0.0	NA